


Fitting your ProSix K9

Your ProSix K9 comes to you assembled. The chest strap (1) has our logo sewn at the front. The rear of the ProSix K9 has two leg straps (2). Please detach these straps; they will hang down centered under the tail when fitting. Now remove the shoulder strap (3) and the hip strap (4) and set them aside placing them so you are aware of where they will go when you add them back to the ProSix K9. The last remaining strap to be removed is the belly strap (5). You are now left with just the chest strap (1) attached to the long body strap (6). This forms a loop that will be placed horizontally on the dog.



MADE IN THE U.S.A.

info@eagleprosix.com

www.eagleprosix.com

The material found on our website and in other promotional materials is provided for educational purposes only and is not to be used for medical advice, diagnosis or treatment. Please seek veterinary advice prior to use as a treatment for rehabilitation or for reconditioning. Neither Eagle Link, LLC nor any authorized retailer of Eagle ProSix will be held liable for any personal accident, injury or damage resulting from the use of the Eagle ProSix. Full responsibility lies with the individual who chooses to use ProSix or ProSix K9 for their horse or dog under his or her care.



ProSix

K9 Fitting Guide

ProSix K9 is the neuromuscular connection that brings the whole dog back online. When your dog first wears ProSix K9, it is important to give the dog time to process the new sensations and the new information between brain and body. A quiet area with as little distraction as possible is best. Please allow the dog to move as it chooses, lying, sitting, rolling, or just standing. Observe the dog for signs of fatigue and allow to rest as needed. We recommend that wear time and activity level be increased gradually. We encourage our clients to take video before using ProSix K9 to have a baseline of movement, posture and symmetry. Video taken over time is great for tracking your dog's progress made from wearing ProSix K9. Reviewing video in slow-motion is optimal for viewing subtle changes in stride, weight distribution and posture.



1. Begin by placing the long body/chest strap around the dog. As seen above, this long body/chest strap loop should be in a straight line from across the chest to across the hind end at the widest part.

Make sure the logo is centered at the front of the chest. The loop should be snug, but not tight. You should be able to slide your hand between the strap and the dog. At this time you can adjust the hoop and loop between the chest and long body straps, one side at a time. Make sure your adjustments overlap equally on both sides.



2. Attach the shoulder strap that you had set aside. (Usually the longer of the two) The shoulder strap attaches on both sides of the dog. Start on one side and lay over the base of the neck and attach to the opposite side.

The shoulder strap should be snug enough to hold the long body strap in a straight line and keep the chest strap in place. Reassess position of straps following activity and adjust as needed.



3. The next step is attaching the hip strap. The hip strap attaches on both sides of the dog. Start on one side and lay over the back of the dog and attach to the opposite side. The hip strap has to be far enough forward that it doesn't slide towards the tail when the dog is in motion. The hip strap keeps the long body strap in a straight line.



4. Now add the belly strap. The belly strap will fit around the dog where the belly starts to get narrower, toward the hind legs. The belly strap should fit over the last true ribs. Begin by attaching the soft loop of the belly band to the long body hook material on one side of the dog. Now go over the dog's back down the opposite side and under the belly.

Continue to bring the belly strap up to the original side and attach it to itself. Check that this strap is snug enough to provide support but never tight. Check that you can easily pull away and return into place (gently, do not snap).



5. Next, find the hind leg straps sewn to the long body strap, under the dog's tail. Pull each strap underneath the hind leg on that side, and bring it up to attach to the long body strap, in front of the hind leg as seen in the picture above. These straps may be crossed to attach on opposite side on female dogs.

Once this process has been completed for your dog and you are happy with the fit, you shouldn't have to go through this entire process again. Simply remove by detaching the hind leg straps and the belly strap. Slide the ProSix K9 up and over the dog's head. To put the ProSix K9 on again, put the dog's head through the body/shoulder strap opening, and pull down over the hind, attach leg straps and belly strap.

